


# em-scalebartools

Fiji macro toolset to quickly add a scale bar with reasonable size to an image. Developed for electron microscopy. | ¶ Please double-check the scale-bar length for possible rounding errors, especially if it shows 1 µm, 1 nm, 1 mm, ... . | |-----|  Example  
Image source: Cropped image of Cobaea scandens pollen by Marie Majaura from Wikipedia.

## Examples


Using QuickScaleBar on a HRTEM image. Note the similar size of the scale bars for the 4096<sup>2</sup> image (center) and the cropped 512<sup>2</sup> ROI image (right).

 Example

Using FEI Crop Scalebar on an SEM image.

 Example


Batch conversion of SEM images ( `Process -> Batch -> Macro...` ) from tiff to png using `FEI_Crop_Scalebar.ijm`.

 Example


## EMScaleBarTools in action

Examples are from v0.2.

Basic usage with cropping of a TFS/FEI databar, addition of a scale bar, moving and removing of the scale bar:

 Example

The next GIF shows the application of `Auto re-scale images` to upscale a small (in pixels) inset of an image:

 Example

The next GIF shows an example workflow when working with presentations (here Microsoft PowerPoint). For a horizontal alignment of images (here an SEM and an HAADF-STEM image) with the same desired image height, the scale bar reference is switched to `Height`. Note the automatic handling of unit-switching and rescaling as in the previous example. The image are then copied via the hotkey c to the system clipboard and pasted into PowerPoint.

 Example

## Macro description

### QuickScaleBar Tool (Icon: **SB**)

- One-click action to add a scale bar to an image. Right click opens the options menu.
- The scale bar height and font size is adjusted based on image height (or width).

- The scale bar width is adjusted based on scaled image width (or height, or larger/smaller of the two) and rounded to next "good looking" number.
- The scale bar appearance can be set up just like the normal settings for `Analyze -> Tools -> Scale bar...`.
- Optional: Automatically switch units to make scale bar more appealing. E.g., an image with horizontal field width of 0.25  $\mu\text{m}$  will be switched to 250 nm. The scale bar will then also be in nm.
- Optional: Automatically re-scale image to (at least) a specified image size in pixels without interpolation (= nearest neighbor interpolation). This is convenient for programs like PowerPoint which like to automatically interpolate "small" images.
- Optional: Run custom macro commands provided in the options menu, e.g. `run('mpl-viridis');` to change LUT to viridis.

## FEI Crop Scalebar Tool (Icon: FEI)

- One-click action to crop away the databar from an FEI/TFS SEM/FIB image and to add a scale bar. Right click opens the options menu.
- Scale bar behaviour is the same as for QuickScaleBar tool and settings are taken from the QuickScaleBar options.
- Optional: Run custom macro commands provided in the options menu, e.g. `run('mpl-viridis');` to change LUT to viridis.
- Especially useful for batch conversion of SEM/FIB images (run from `Process -> Batch -> Macro...`): In the batch processing menu insert the macro command `runMacro('FEI_Crop_Scalebar.ijm');`.

For a description of the other tools, take a look at the [wiki](#). 

## Requirements and Installation

- Cropping the FEI/TFS info bar requires the useful [EM tool](#) plugin by **IMBalENce** as FEI/TFS images are scaled with [SEM FEI metadata scale](#). Install via the Fiji update site.
- Download the latest [release](#), extract the `macros` folder, and copy it to your Fiji installation folder. It will add the `FEI_Crop_Scalebar.ijm` macro to the macros folder and the `EMScaleBarTools.ijm` toolset to the `macros/toolset` folder.
- Restart Fiji and select the `EMScaleBarTools` from `More Tools...` (>>) menu.
- Legacy option for v0.2 and below: [Adding new FEI/TFS microscope types](#)

## Documentation

# **Changelog**

## **Other useful scalebar tools**

- Python: [matplotlib-scalebar](#) by ppinard
- DM/GMS: [Scale Bar Control](#) by D. R. G. Mitchell
- Fiji/ImageJ: [asc-ImageJ-Fancy-Labels](#) by peterjlee
- Fiji/ImageJ: [Scale Bar Tools for Microscopes](#) by Gilles Carpentier

## **Citing**

If you want, you can cite this project via Zenodo:

DOI [10.5281/zenodo.7799311](https://doi.org/10.5281/zenodo.7799311)